

For lab use only.

Cast Date: _____

Doctors name: _____ Patients last name: _____ First: _____

Address: _____ Age: _____ Sex: _____ Weight: _____ Shoe size: _____

City: _____ Province: _____ Phone: _____ Fax: _____

Diagnosis: _____

Orthotic Selection

- Shock
- Soft
- Youth
- ucbl
- Sport Flex
- All Sport
- Casual Ultra Flex
- Active Diabetic
- Dress: Style Pump Flex

Shell Material Selection

- White Polypropylene 2 mm
- White Polypropylene 3 mm
- Carbon Flex 2mm

Posting Instructions

Forefoot

- Lab Discretion
- FF Intrinsic
- FF Extrinsic
- Varus Left _____ Right _____
- Valgus Left _____ Right _____

Rearfoot

- RF Intrinsic in Shell
- RF Extrinsic
- Heel Lift Left _____ mm Right _____ mm
- Varus Left _____ Right _____
- Valgus Left _____ Right _____

Castwork

- Low Arch Fill
- Standard Arch Fill
- No Arch Fill

Heel Cup Depth

- Deep
- Shallow

Top Cover Length

- Meta Length
- Sulcus Length
- Full Length

Top Cover Material

- Blue ETC 1/8" Pad
- Black ETC 1/8" Pad
- Vinyl 1/16" Pad 1/8" Pad
- Star Suede
- Suede
- Leather

Accommodations

- | | R | L |
|--------------------|--------------------------|--------------------------|
| First Ray Cut-out | <input type="checkbox"/> | <input type="checkbox"/> |
| Arch Pad | <input type="checkbox"/> | <input type="checkbox"/> |
| Heel Pad | <input type="checkbox"/> | <input type="checkbox"/> |
| Heel Spur Pad | <input type="checkbox"/> | <input type="checkbox"/> |
| Meta Bar | <input type="checkbox"/> | <input type="checkbox"/> |
| Meta Pad | <input type="checkbox"/> | <input type="checkbox"/> |
| Lateral Flange | <input type="checkbox"/> | <input type="checkbox"/> |
| Medial Flange | <input type="checkbox"/> | <input type="checkbox"/> |
| Arch Fill - Soft | <input type="checkbox"/> | <input type="checkbox"/> |
| Arch Fill - Firm | <input type="checkbox"/> | <input type="checkbox"/> |
| Morton's Extension | <input type="checkbox"/> | <input type="checkbox"/> |
| Neuroma Pad | <input type="checkbox"/> | <input type="checkbox"/> |

Special Instructions

Mark Accommodations on cast and below

Office Supplies

- Order Forms
- Shipping Boxes

Bio-Foam

- 6 Pack 44 + \$8.00 shipping
- 25 Pack 44 + \$20.00 shipping
- 50 Pack 44 + \$20.00 shipping

Special Instructions:



Right

Left

Shoe Selection

Shoe : _____

Size : _____ Colour : _____

Width : _____ with orthotic shoe only

Shoe : _____

Size : _____ Colour : _____

Width : _____ with orthotic shoe only

Orthotics

Shock

Device designed to provide high levels of control with design features allowing maximum shock absorption. Best suited for casual, workplace and active footwear which allow depth for device.

- High control level
- Semi rigid poly pro shell with arch fill
- Full length puff top covering with nyplex bottom cover

Soft

Trilaminar shell, designed to provide soft cushioning and enhanced support. Ideal device for sensitive clients who require gentle support. Full fitting casual footwear suggested. Not recommended for dress or athletic footwear due to bulk and design.

Youth

Orthotics constructed utilizing a semi-flexible shell design to provide motion control and comfort for children. Suited for athletic footwear and casual day to day youth footwear.

- 2mm carbon flex or poly pro shell design
- Sport ETC top covering

ucbl

Specialty device for children and adults with flexible flat feet requiring maximum control and correction. Rigid device with medial and lateral phlanges and comfort eva top covering.

- RCH or polypro shell design
- Raised medial and lateral phlanges
- Puff top covering

Sport Flex

Designed to meet the needs of active clients requiring motion control, and shock absorption. This is an ideal device for athletic activities and workplace footwear for individuals who perform long periods of standing or walking..

- 1mm or 2mm carbon shell
- Poron arch fill
- Sport ETC top cover
- Nyplex bottom covering

All Sport

Specifically constructed to meet the needs of patients who require increased motion control in a lightweight resilient device. Fits most athletic footwear types.

- 2mm or 3mm poly pro shell
- Full length sport etc coverings

Casual Ultra Flex

Orthotics designed for patients who require support and flexibility in their day to day activities. Fits most shoes with the exception of narrow dress style footwear. An excellent product for clients who require one pair of orthotics for use in both athletic or casual footwear.

- 2mm or 3mm carbon flex shell
- Full length supreme etc covering

Dress - Style

For Him, or Her, low profile device for use in slip on or lace dress shoes, designed to promote enhanced biomechanical control with reduced bulk.

- 1mm or 2mm carbon shell
- 1/16 starsuede top covering
- Narrow orthotic grind

Dress - Pump Flex

Sleek low profile orthotic allowing mild biomechanical control for use in women's high heel shoes.

- 1mm carbon shell
- Narrow tapered grind design
- 1/16 starsuede top covering

Active Diabetic

For the diabetic patient who's primary fitness activity is walking. Ideal device for the type 2 patient who requires stability, support and palliation during their daily activities.

- 1mm or 2mm carbon flex shell
- Poron arch fill
- Poron 1/8 and plastizote top covering